

IHE Integration Statement		Date
Vendor	Product Name	Version
NIHON KOHDEN	PrimeVita/PrimeVitaPlus	Ver.02-11/Ver.02-03
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:		
Integration Profiles Implemented	Actors Implemented	Options Implemented
Radiology		
Scheduled Workflow	Order Filler	Japanese
	PPS Manager	Japanese
Patient Information Reconciliation	Order Filler	Japanese
	PPS Manager	Japanese
Cardiology		
ECHO	Order Filler	Japanese
	Image Manager/Archive	Japanese
	Image Display	Japanese
IVI	Order Filler	Japanese
	PPS Manager	Japanese
	Image Manager/Archive	Japanese
	Image Display	Japanese
CATH	Order Filler	Japanese
	Image Manager/Archive	Japanese
	Image Display	Japanese
STRESS	Order Filler	Japanese
	Image Manager/Archive	Japanese
	Image Display	Japanese
ECG	Display	Japanese
	Information Source	Japanese
ITI		
CT	Time Client	
PAM	Patient Demographics Consumer	
	Patient Encounter Consumer	
Endoscopy		
EWF	Endoscopy Order Filler	
	Execution Information Creator	
Internet address for vendor's IHE information:		
http://www.nihonkohden.co.jp/iryu/products/examroomsys/01_dis/prm3000.html		
http://www.nihonkohden.co.jp/iryu/products/examroomsys/01_dis/prm4000.html		
Links to Standards Conformance Statements for the Implementation		
HL7	Not Applicable	
DICOM	http://www.nihonkohden.co.jp/iryu/products/examroomsys/01_dis/prm3000.html	
	http://www.nihonkohden.co.jp/iryu/products/examroomsys/01_dis/prm4000.html	
Links to general information on IHE		
In North America:	In Europe:	In Japan:
http://www.ihe.net/	http://www.ihe-europe.org/	http://www.ihe-j.org/

- *IHE-J 2008 Connect-a-thon 2008/10/27-31
- *IHE-J 2009 Connect-a-thon 2009/10/26-30
- *IHE-J 2010 Connect-a-thon 2010/10/18-22
- *IHE-J 2011 Connect-a-thon 2011/11/7-11
- *IHE-J 2012 Connect-a-thon 2012/10/29-11/2
- *IHE-J 2013 Connect-a-thon 2013/10/28-11/1